

— ABLYBIO, Help Your Research



# HNF-4-alpha (YD13834) Rabbit mAb

货号: **AYD13230**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	53kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Nucleus
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a nuclear transcription factor which binds DNA as a homodimer. The encoded protein controls the expression of several genes, including hepatocyte nuclear factor 1 alpha, a transcription factor which regulates the expression of several hepatic genes. This gene may play a role in development of the liver, kidney, and intestines. Mutations in this gene have been associated with monogenic autosomal dominant non-insulin-dependent diabetes mellitus type I. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms.
基因ID	15378
基因名	Hnf4a
Swiss	P49698 ( <a href="https://www.uniprot.org/uniprotkb/P49698/entry">https://www.uniprot.org/uniprotkb/P49698/entry</a> )
别名	HNF-4-alpha (YD13834),HNF-4-alpha (YD13834) Rabbit mAb,Hnf4a,Nuclear receptor subfamily 2 group A member 1,Transcription factor 14,Transcription factor HNF-4,Hnf-4,Hnf4,Nr2a1,Tcf14

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)