

Xanthine Oxidase (YD13832) Rabbit mAb

货号: **AYD13232**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	146kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,Mouse liver,Rat lung,Rat liver
细胞定位	Cytoplasm, Peroxisome, Secreted
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	Xanthine dehydrogenase belongs to the group of molybdenum-containing hydroxylases involved in the oxidative metabolism of purines. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. Xanthine dehydrogenase can be converted to xanthine oxidase by reversible sulfhydryl oxidation or by irreversible proteolytic modification. Defects in xanthine dehydrogenase cause xanthinuria, may contribute to adult respiratory stress syndrome, and may potentiate influenza infection through an oxygen metabolite-dependent mechanism.
基因ID	7498
基因名	XDH

Swiss	P47989
别名	Xanthine Oxidase (YD13832)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn