

— ABLYBIO, Help Your Research



HDAC6 (YD33228) Rabbit mAb

货号: **AYD13243**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	131kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat liver
细胞定位	Cytoplasm, cytoskeleton, Nucleus, Perikaryon, Cell projection, dendrite, axon, cilium, microtubule organizing center, centrosome, cilium basal body
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to class II of the histone deacetylase/acuc/apha family. It contains an internal duplication of two catalytic domains which appear to function independently of each other. This protein possesses histone deacetylase activity and represses transcription.
基因ID	10013
基因名	HDAC6
Swiss	Q9UBN7 (https://www.uniprot.org/uniprotkb/Q9UBN7/entry)
别名	HDAC6 (YD33228),HDAC6 (YD33228) Rabbit mAb,HDAC6,E3 ubiquitin-protein ligase HDAC6,Tubulin-lysine deacetylase HDAC6,KIAA0901

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)