

MYBBP1A (YD33216) Rabbit mAb

货号: **AYD13245**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	149kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SW620,MCF7,Mouse liver,Rat liver
细胞定位	Cytoplasm, Nucleus, nucleolus
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a nucleolar transcriptional regulator that was first identified by its ability to bind specifically to the Myb proto-oncogene protein. The encoded protein is thought to play a role in many cellular processes including response to nucleolar stress, tumor suppression and synthesis of ribosomal DNA. Alternate splicing results in multiple transcript variants.
基因ID	10514
基因名	MYBBP1A
Swiss	Q9BQG0

别名	MYBBP1A (YD33216)
----	-------------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn