

— ABLYBIO, Help Your Research



# HOPX (YD14257) Rabbit mAb

货号: **AYD13250**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC-P:</b> 1:100 - 1:200 <b>IP:</b> 1:20 - 1:50
理论分子量	8kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus, Cytoplasm
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Atypical homeodomain protein which does not bind DNA and is required to modulate cardiac growth and development. Acts via its interaction with SRF, thereby modulating the expression of SRF-dependent cardiac-specific genes and cardiac development. Prevents SRF-dependent transcription either by inhibiting SRF binding to DNA or by recruiting histone deacetylase (HDAC) proteins that prevent transcription by SRF. Overexpression causes cardiac hypertrophy (PubMed:12297045, PubMed:12297046). Acts as a co-chaperone for HSPA1A and HSPA1B chaperone proteins and assists in chaperone-mediated protein refolding (By similarity)
基因ID	3134
基因名	Hopx
Swiss	Q8R1H0 ( <a href="https://www.uniprot.org/uniprotkb/Q8R1H0/entry">https://www.uniprot.org/uniprotkb/Q8R1H0/entry</a> )
别名	HOPX (YD14257),HOPX (YD14257) Rabbit mAb,Hopx,Homeobox-only protein,Odd homeobox protein 1,Hod,Hop,Ob1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)