

— ABLYBIO, Help Your Research



# Phospho-Rb (Ser780) (YD14539) Rabbit mAb

货号: **AYD13253**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	106kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-431
细胞定位	Nucleus, Cytoplasm
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a negative regulator of the cell cycle and was the first tumor suppressor gene found. The encoded protein also stabilizes constitutive heterochromatin to maintain the overall chromatin structure. The active, hypophosphorylated form of the protein binds transcription factor E2F1. Defects in this gene are a cause of childhood cancer retinoblastoma (RB), bladder cancer, and osteogenic sarcoma.
基因ID	5925
基因名	RB1
Swiss	P06400 ( <a href="https://www.uniprot.org/uniprotkb/P06400/entry">https://www.uniprot.org/uniprotkb/P06400/entry</a> )
别名	Phospho-Rb (Ser780) (YD14539), Phospho-Rb (Ser780) (YD14539) Rabbit mAb, RB1, p105-Rb, p110-RB1, pRb, pp110

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)