

— ABLYBIO, Help Your Research



# Kv4.2/KCND2 (YD33065) Rabbit mAb

货号: **AYD13266**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC FC
推荐浓度	
理论分子量	71kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	Cell membrane, Cell projection, dendrite, Synapse, Perikaryon, Postsynaptic cell membrane, dendritic spine, Cell junction
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Voltage-gated potassium (Kv) channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. Four sequence-related potassium channel genes - shaker, shaw, shab, and shal - have been identified in Drosophila, and each has been shown to have human homolog(s). This gene encodes a member of the potassium channel, voltage-gated, shal-related subfamily, members of which form voltage-activated A-type potassium ion channels and are prominent in the repolarization phase of the action potential. This member mediates a rapidly inactivating, A-type outward potassium current which is not under the control of the N terminus as it is in Shaker channels.
基因ID	3751
基因名	KCND2
Swiss	Q9NZV8 ( <a href="https://www.uniprot.org/uniprotkb/Q9NZV8/entry">https://www.uniprot.org/uniprotkb/Q9NZV8/entry</a> )
别名	Kv4.2/KCND2 (YD33065),Kv4.2/KCND2 (YD33065) Rabbit mAb,KCND2,Potassium voltage-gated channel subfamily D member 2,Voltage-gated potassium channel subunit Kv4.2,KIAA1044

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)