

— ABLYBIO, Help Your Research



CD11b (YD33044) Rabbit mAb

货号: **AYD13268**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	127kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spleen,C2C12(Negative control),RAW264.7
细胞定位	Cell membrane, Membrane raft
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	Enables complement component C3b binding activity; heparan sulfate proteoglycan binding activity; and heparin binding activity. Contributes to cargo receptor activity. Involved in several processes, including central nervous system development; endocytosis; and positive regulation of neuron death. Acts upstream of or within several processes, including activated T cell proliferation; leukocyte migration; and microglia development. Located in external side of plasma membrane and nucleus. Is integral component of membrane. Part of integrin alphaM-beta2 complex. Is expressed in several structures, including central nervous system; heart; lymphatic vessel; thymus; and yolk sac. Human ortholog(s) of this gene implicated in lupus nephritis. Orthologous to human ITGAM (integrin subunit alpha M).
基因ID	16409
基因名	Itgam
Swiss	P05555 (https://www.uniprot.org/uniprotkb/P05555/entry)
别名	CD11b (YD33044),CD11b (YD33044) Rabbit mAb,Itgam,CD11 antigen-like family member B,CR-3 alpha chain,Cell surface glycoprotein MAC-1 subunit alpha,Leukocyte adhesion receptor MO1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)