

— ABLYBIO, Help Your Research



PON1 (YD33003) Rabbit mAb

货号: **AYD13286**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	40kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Jurkat,Mouse thymus,Mouse heart,Rat liver
细胞定位	Secreted, extracellular space
纯化	亲和纯化

抗原信息

抗原信息	
------	--


靶点信息

研究背景	The enzyme encoded by this gene is an arylesterase that mainly hydrolyzes paroxon to produce p-nitrophenol. Paroxon is an organophosphorus anticholinesterase compound that is produced in vivo by oxidation of the insecticide parathion. Polymorphisms in this gene are a risk factor in coronary artery disease. The gene is found in a cluster of three related paraoxonase genes at 7q21.3.
------	--

基因ID	5444
基因名	PON1
Swiss	P27169 (https://www.uniprot.org/uniprotkb/P27169/entry)
别名	PON1 (YD33003),PON1 (YD33003) Rabbit mAb,PON1,Aromatic esterase 1,K-45,Serum aryldialkylphosphatase 1,PON

产品验证

Human plasma



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)