

— ABLYBIO, Help Your Research



# Ctip2 (YD32960) Rabbit mAb

货号: **AYD13312**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	96kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse thymus
细胞定位	Nucleus
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a C2H2-type zinc finger protein and is closely related to BCL11A, a gene whose translocation may be associated with B-cell malignancies. Although the specific function of this gene has not been determined, the encoded protein is known to be a transcriptional repressor, and is regulated by the NURD nucleosome remodeling and histone deacetylase complex. Four alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
基因ID	64919
基因名	BCL11B
Swiss	Q9C0K0 ( <a href="https://www.uniprot.org/uniprotkb/Q9C0K0/entry">https://www.uniprot.org/uniprotkb/Q9C0K0/entry</a> )
别名	Ctip2 (YD32960),Ctip2 (YD32960) Rabbit mAb,BCL11B,B-cell CLL/lymphoma 11B,COUP-TF-interacting protein 2,Radiation-induced tumor suppressor gene 1 protein,CTIP2,RIT1

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)