

— ABLYBIO, Help Your Research



Cytokeratin 6 (YD32690) Rabbit mAb

货号: **AYD13381**

产品信息

| | |
|-------|---|
| 反应 | Human |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC-P ICC/IF |
| 推荐浓度 | |
| 理论分子量 | 60kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | A-431,293T |
| 细胞定位 | |
| 纯化 | 亲和纯化 |



抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|---|
| 研究背景 | The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. As many as six of this type II cytokeratin (KRT6) have been identified; the multiplicity of the genes is attributed to successive gene duplication events. The genes are expressed with family members KRT16 and/or KRT17 in the filiform papillae of the tongue, the stratified epithelial lining of oral mucosa and esophagus, the outer root sheath of hair follicles, and the glandular epithelia. This KRT6 gene in particular encodes the most abundant isoform. Mutations in these genes have been associated with pachyonychia congenita. In addition, peptides from the C-terminal region of the protein have antimicrobial activity against bacterial pathogens. The type II cytokeratins are clustered in a region of chromosome 12q12-q13. |
| 基因ID | 3853 |
| 基因名 | KRT6A |
| Swiss | P02538 (https://www.uniprot.org/uniprotkb/P02538/entry) |
| 别名 | Cytokeratin 6 (YD32690),Cytokeratin 6 (YD32690) Rabbit mAb,KRT6A,Cytokeratin-6A,Cytokeratin-6D,Keratin-6A,Type-II keratin Kb6,K6A,KRT6D |

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)