

— ABLYBIO, Help Your Research



Cpn10 (YD12474) Rabbit mAb

货号: **AYD13424**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	11kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Mouse liver,Mouse kidney,Mouse heart,Rat liver
细胞定位	Mitochondrion matrix
纯化	亲和纯化




抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a major heat shock protein which functions as a chaperonin. Its structure consists of a heptameric ring which binds to another heat shock protein in order to form a symmetric, functional heterodimer which enhances protein folding in an ATP-dependent manner. This gene and its co-chaperonin, HS PD1, are arranged in a head-to-head orientation on chromosome 2. Naturally occurring read-through transcription occurs between this locus and the neighboring locus MOBKL3.
基因ID	3336
基因名	HSPE1
Swiss	P61604 (https://www.uniprot.org/uniprotkb/P61604/entry)
别名	Cpn10 (YD12474),Cpn10 (YD12474) Rabbit mAb,HSPE1,10 kDa chaperonin,Chaperonin 10,Early-pregnancy factor,Heat shock protein family E member 1

产品验证





实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)