

HIF1AN (YD31952) Rabbit mAb

货号: **AYD13431**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	40kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,K-562,U-937,293T,Mouse heart,Mouse skeletal muscle
细胞定位	Nucleus, Cytoplasm, perinuclear region
纯化	

抗原信息





抗原信息	
------	--

靶点信息

研究背景	Hydroxylates HIF-1 alpha at 'Asn-803' in the C-terminal transactivation domain (CAD. Functions as an oxygen sensor and, under normoxic conditions, the hydroxylation prevents interaction of HIF-1 with transcriptional coactivators including Cbp/p300-interacting transactivator. Involved in transcriptional repression through interaction with HIF1A, VHL and histone deacetylases. Hydroxylates specific Asn residues within ankyrin repeat domains (ARD of NFKB1, NFKBIA, NOTCH1, ASB4, PPP1R12A and several other ARD-containing proteins. Also hydroxylates Asp and His residues within ARDs of ANK1 and TNKS2, respectively. Negatively regulates NOTCH1 activity, accelerating myogenic differentiation. Positively regulates ASB4 activity, promoting vascular differentiation.
基因ID	55662

基因名	HIF1AN
Swiss	Q9NWT6
别名	HIF1AN (YD31952)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn