

— ABLYBIO, Help Your Research



# Cytokeratin 4 (YD11986) Rabbit mAb

货号: **AYD13441**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	56kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,A-431,Mouse stomach
细胞定位	
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in differentiated layers of the mucosal and esophageal epithelia with family member KRT13. Mutations in these genes have been associated with White Sponge Nevus, characterized by oral, esophageal, and anal leukoplakia. The type II cytokeratins are clustered in a region of chromosome 12q12-q13. [provided by RefSeq, Jul 2008]
基因ID	3851
基因名	KRT4
Swiss	P19013 ( <a href="https://www.uniprot.org/uniprotkb/P19013/entry">https://www.uniprot.org/uniprotkb/P19013/entry</a> )
别名	Cytokeratin 4 (YD11986),Cytokeratin 4 (YD11986) Rabbit mAb,KRT4,Cytokeratin-4,Keratin-4,Type-II keratin Kb4,CYK4

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)