

— ABLYBIO, Help Your Research



# SERPING1 (YD32550) Rabbit mAb

货号: **AYD13444**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	55kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human plasma
细胞定位	Secreted
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--


## 靶点信息

研究背景	This gene encodes a highly glycosylated plasma protein involved in the regulation of the complement cascade. Its protein inhibits activated C1r and C1s of the first complement component and thus regulates complement activation. Deficiency of this protein is associated with hereditary angioneurotic oedema (HANE). Alternative splicing results in multiple transcript variants encoding the same isoform.
------	---

基因ID	710
基因名	SERPING1
Swiss	P05155 ( <a href="https://www.uniprot.org/uniprotkb/P05155/entry">https://www.uniprot.org/uniprotkb/P05155/entry</a> )
别名	SERPING1 (YD32550),SERPING1 (YD32550) Rabbit mAb,SERPING1,C1 esterase inhibitor,C1-inhibiting factor,Serpin G1,C1IN,C1NH

## 产品验证

Human plasma



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)