

CCL27 (YD32526) Rabbit mAb

货号: **AYD13455**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	13kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted, Nucleus
纯化	


抗原信息

抗原信息	
------	--

靶点信息

研究背景	Chemotactic factor that attracts skin-associated memory T-lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. May play a role in cell migration during embryogenesis. Nuclear forms may facilitate cellular migration by inducing cytoskeletal relaxation. Binds to CCR10
基因ID	3134
基因名	Ccl27
Swiss	Q9Z1X0
别名	CCL27 (YD32526)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn