

# PIM1 (YD14753) Rabbit mAb

货号: **AYD13460**

## 产品信息

|       |   |
|-------|---|
| 反应    | Human,Mouse,Rat   |
| 宿主    | Rabbit  |
| 克隆性   | Monoclonal  |
| 预测反应  |   |
| 应用    | WB IHC-P ICC/IF FC  |
| 推荐浓度  |   |
| 理论分子量 | 36kDa   |
| 实测分子量 |   |
| 形式    | Liquid  |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated  |
| 阳性对照  | H460,Mouse kidney,Rat kidney  |
| 细胞定位  | Cytoplasm, Nucleus, Cell membrane   |
| 纯化    |   |

## 抗原信息



|      |  |
|------|--|
| 抗原信息 |  |
|------|--|

## 靶点信息

|      |   |
|------|---|
| 研究背景 | The protein encoded by this gene belongs to the Ser/Thr protein kinase family, and PIM subfamily. This gene is expressed primarily in B-lymphoid and myeloid cell lines, and is overexpressed in hematopoietic malignancies and in prostate cancer. It plays a role in signal transduction in blood cells, contributing to both cell proliferation and survival, and thus provides a selective advantage in tumorigenesis. Both the human and orthologous mouse genes have been reported to encode two isoforms (with preferential cellular localization) resulting from the use of alternative in-frame translation initiation codons, the upstream non-AUG (CUG) and downstream AUG codons (PMIDs:16186805, 1825810). |
| 基因ID | 5292  |
| 基因名  | PIM1  |

|       |                |
|-------|----------------|
| Swiss | P11309         |
| 别名    | PIM1 (YD14753) |

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)