

— ABLYBIO, Help Your Research



# Cytokeratin 14 (YD32496) Rabbit mAb

货号: **AYD13465**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse brain,Rat brain
细胞定位	Cytoplasm, Nucleus
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the keratin family, the most diverse group of intermediate filaments. This gene product, a type I keratin, is usually found as a heterotetramer with two keratin 5 molecules, a type II keratin. Together they form the cytoskeleton of epithelial cells. Mutations in the genes for these keratins are associated with epidermolysis bullosa simplex. At least one pseudogene has been identified at 17p12-p11.
基因ID	3861
基因名	KRT14
Swiss	P02533 ( <a href="https://www.uniprot.org/uniprotkb/P02533/entry">https://www.uniprot.org/uniprotkb/P02533/entry</a> )
别名	Cytokeratin 14 (YD32496),Cytokeratin 14 (YD32496) Rabbit mAb,KRT14,Cytokeratin-14,Keratin-14

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)