

— ABLYBIO, Help Your Research



# CD276 (YD11447) Rabbit mAb

货号: **AYD13477**

## 产品信息

反应	Human, Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC IP
推荐浓度	
理论分子量	34kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Molt-4,Jurkat,Mouse spleen,Mouse thymus,Rat lung
细胞定位	Membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene belongs to the immunoglobulin superfamily, and thought to participate in the regulation of T-cell-mediated immune response. Studies show that while the transcript of this gene is ubiquitously expressed in normal tissues and solid tumors, the protein is preferentially expressed only in tumor tissues. Additionally, it was observed that the 3' UTR of this transcript contains a target site for miR29 microRNA, and there is an inverse correlation between the expression of this protein and miR29 levels, suggesting regulation of expression of this gene product by miR29. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	80381
基因名	Cd276
Swiss	Q8VE98 ( <a href="https://www.uniprot.org/uniprotkb/Q8VE98/entry">https://www.uniprot.org/uniprotkb/Q8VE98/entry</a> )
别名	CD276 (YD11447),CD276 (YD11447) Rabbit mAb,Cd276,B7 homolog 3,Costimulatory molecule,B7h3

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)