

— ABLYBIO, Help Your Research



# Periostin (YD14125) Rabbit mAb

货号: **AYD13490**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	93kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Golgi apparatus, Secreted, extracellular space, extracellular matrix
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a secreted extracellular matrix protein that functions in tissue development and regeneration, including wound healing, and ventricular remodeling following myocardial infarction. The encoded protein binds to integrins to support adhesion and migration of epithelial cells. This protein plays a role in cancer stem cell maintenance and metastasis. Mice lacking this gene exhibit cardiac valve disease, and skeletal and dental defects. Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	10631
基因名	POSTN
Swiss	Q15063 ( <a href="https://www.uniprot.org/uniprotkb/Q15063/entry">https://www.uniprot.org/uniprotkb/Q15063/entry</a> )
别名	Periostin (YD14125), Periostin (YD14125) Rabbit mAb, POSTN, Osteoblast-specific factor 2, OSF2

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)