

— ABLYBIO, Help Your Research



# UFM1 (YD11381) Rabbit mAb

货号: **AYD13504**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	9kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat pancreas
细胞定位	Nucleus, Cytoplasm
纯化	亲和纯化

## 抗原信息




抗原信息	
------	--

## 靶点信息

研究背景	UFM1 is a ubiquitin-like protein that is conjugated to target proteins by E1-like activating enzyme UBA5 (UBE1DC1; MIM 610552) and E2-like conjugating enzyme UFC1 (MIM 610554) in a manner analogous to ubiquitylation (see UBE2M; MIM 603173) (Komatsu et al., 2004 [PubMed 15071506]).
------	---

基因ID	51569
基因名	UFM1
Swiss	P61960 ( <a href="https://www.uniprot.org/uniprotkb/P61960/entry">https://www.uniprot.org/uniprotkb/P61960/entry</a> )
别名	UFM1 (YD11381),UFM1 (YD11381) Rabbit mAb,UFM1,C13orf20

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)