

— ABLYBIO, Help Your Research



p27 KIP 1 (YD32366) Rabbit mAb

货号: **AYD13505**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	22kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus, Cytoplasm, Endosome
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Important regulator of cell cycle progression (PubMed:12972555, PubMed:8033213). Inhibits the kinase activity of CDK2 bound to cyclin A, but has little inhibitory activity on CDK2 bound to SPDYA (By similarity). Involved in G1 arrest. Potent inhibitor of cyclin E- and cyclin A-CDK2 complexes (PubMed:8033213). Forms a complex with cyclin type D-CDK4 complexes and is involved in the assembly, stability, and modulation of CCND1-CDK4 complex activation. Acts either as an inhibitor or an activator of cyclin type D-CDK4 complexes depending on its phosphorylation state and/or stoichiometry
基因ID	3134
基因名	Cdkn1b
Swiss	P46414 (https://www.uniprot.org/uniprotkb/P46414/entry)
别名	p27 KIP 1 (YD32366), p27 KIP 1 (YD32366) Rabbit mAb, Cdkn1b, Cyclin-dependent kinase inhibitor p27, p27 Kip1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)