

— ABLYBIO, Help Your Research



Syntaxin Binding Protein 1 (YD33593) Rabbit mAb

货号: AYD13510

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	68kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,K-562,Mouse lung,Mouse brain
细胞定位	Cytoplasm, cytosol, Membrane
纯化	亲和纯化




抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a syntaxin-binding protein. The encoded protein appears to play a role in release of neurotransmitters via regulation of syntaxin, a transmembrane attachment protein receptor. Mutations in this gene have been associated with infantile epileptic encephalopathy-4. Alternatively spliced transcript variants have been described.
基因ID	6812
基因名	STXBP1
Swiss	P61764 (https://www.uniprot.org/uniprotkb/P61764/entry)
别名	Syntaxin Binding Protein 1 (YD33593),Syntaxin Binding Protein 1 (YD33593) Rabbit mAb,STXBP1,MUNC18-1,N-Sec1,Protein unc-18 homolog 1,Protein unc-18 homolog A,p67,UNC18A

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)