

— ABLYBIO, Help Your Research



PU.1/Spi1 (YD14729) Rabbit mAb

货号: **AYD13554**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	31kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-937,Mouse spleen
细胞定位	Nucleus
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes an ETS-domain transcription factor that activates gene expression during myeloid and B-lymphoid cell development. The nuclear protein binds to a purine-rich sequence known as the PU-box found near the promoters of target genes, and regulates their expression in coordination with other transcription factors and cofactors. The protein can also regulate alternative splicing of target genes. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	6688
基因名	SPI1
Swiss	P17947 (https://www.uniprot.org/uniprotkb/P17947/entry)
别名	PU.1/Spi1 (YD14729),PU.1/Spi1 (YD14729) Rabbit mAb,SPI1,31 kDa-transforming protein

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)