

— ABLYBIO, Help Your Research



# CBX4 (YD18870) Rabbit mAb

货号: **AYD13563**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	61kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,HepG2,SW620,Jurkat,Mouse brain,Mouse liver,Rat liver
细胞定位	Nucleus, Nucleus speckle
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	E3 SUMO-protein ligase which facilitates SUMO1 conjugation by UBE2I. Involved in the sumoylation of HN RNPk, a p53/TP53 transcriptional coactivator, hence indirectly regulates p53/TP53 transcriptional activation resulting in p21/CDKN1A expression. Monosumoylates ZNF131.
------	--

基因ID	8535
基因名	CBX4
Swiss	O00257 ( <a href="https://www.uniprot.org/uniprotkb/O00257/entry">https://www.uniprot.org/uniprotkb/O00257/entry</a> )
别名	CBX4 (YD18870),CBX4 (YD18870) Rabbit mAb,CBX4,Chromobox protein homolog 4,Polycomb 2 homolog

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)