

LCAT (YD32142) Rabbit mAb

货号: **AYD13564**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human plasma
细胞定位	Secreted
纯化	


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes the extracellular cholesterol esterifying enzyme, lecithin-cholesterol acyltransferase. The esterification of cholesterol is required for cholesterol transport. Mutations in this gene have been found to cause fish-eye disease as well as LCAT deficiency.
基因ID	3931
基因名	LCAT
Swiss	P04180
别名	LCAT (YD32142)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn