

— ABLYBIO, Help Your Research



Skp1 (YD32122) Rabbit mAb

货号: **AYD13572**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC
推荐浓度	
理论分子量	19kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,A-431,HeLa,LO2,SW480,K-562,Mouse kidney
细胞定位	
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a component of SCF complexes, which are composed of this protein, cullin 1, a ring-box protein, and one member of the F-box family of proteins. This protein binds directly to the F-box motif found in F-box proteins. SCF complexes are involved in the regulated ubiquitination of specific protein substrates, which targets them for degradation by the proteasome. Specific F-box proteins recognize different target protein(s), and many specific SCF substrates have been identified including regulators of cell cycle progression and development. Studies have also characterized the protein as an RNA polymerase II elongation factor. Alternative splicing of this gene results in two transcript variants. A related pseudogene has been identified on chromosome 7.
基因ID	6500
基因名	SKP1
Swiss	P63208 (https://www.uniprot.org/uniprotkb/P63208/entry)
别名	Skp1 (YD32122),Skp1 (YD32122) Rabbit mAb,SKP1,Cyclin-A/CDK2-associated protein p19,Organ of Corti protein 2,Organ of Corti protein II,RNA polymerase II elongation factor-like protein,SIII,Transcription elongation factor B polypeptide 1-like,p19skp1,EMC19,OCP2,SKP1A

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)