

— ABLYBIO, Help Your Research



# Chromogranin A (YD19749) Rabbit mAb

货号: **AYD13580**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	51kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,mouse kidney
细胞定位	Secreted, Cytoplasmic vesicle, secretory vesicle, neuronal dense core vesicle
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. It is found in secretory vesicles of neurons and endocrine cells. This gene product is a precursor to three biologically active peptides; vasostatin, pancreastatin, and parastatin. These peptides act as autocrine or paracrine negative modulators of the neuroendocrine system. Two other peptides, catestatin and chromofungin, have antimicrobial activity and antifungal activity, respectively. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	1113
基因名	CHGA
Swiss	P10645 ( <a href="https://www.uniprot.org/uniprotkb/P10645/entry">https://www.uniprot.org/uniprotkb/P10645/entry</a> )
别名	Chromogranin A (YD19749),Chromogranin A (YD19749) Rabbit mAb,CHGA,Pituitary secretory protein I,Vasostatin I,Vasostatin II,SL21

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))