

— ABLYBIO, Help Your Research



Septin 8 (YD31950) Rabbit mAb

货号: **AYD13583**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	56kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,HeLa,PC-3(negative),Mouse brain,Mouse heart,Rat brain,Rat lung
细胞定位	Cytoplasm, cytoskeleton, Synapse, Cell projection, axon, Cytoplasmic vesicle, secretory vesicle, synaptic vesicle membrane, Presynapse
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene is a member of the septin family of nucleotide binding proteins, originally described in yeast as cell division cycle regulatory proteins. Septins are highly conserved in yeast, Drosophila, and mouse, and appear to regulate cytoskeletal organization. Disruption of septin function disturbs cytokinesis and results in large multinucleate or polyploid cells. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	23176
基因名	SEPTIN8
Swiss	Q92599 (https://www.uniprot.org/uniprotkb/Q92599/entry)
别名	Septin 8 (YD31950),Septin 8 (YD31950) Rabbit mAb,SEPTIN8,KIAA0202,SEPT8

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)