

Septin 8 (YD31950) Rabbit mAb

货号: **AYD13583**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	56kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,HeLa,PC-3(negative),Mouse brain,Mouse heart,Rat brain,Rat lung
细胞定位	Cytoplasm, cytoskeleton, Synapse, Cell projection, axon, Cytoplasmic vesicle, secretory vesicle, synaptic vesicle membrane, Presynapse
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	This gene is a member of the septin family of nucleotide binding proteins, originally described in yeast as cell division cycle regulatory proteins. Septins are highly conserved in yeast, Drosophila, and mouse, and appear to regulate cytoskeletal organization. Disruption of septin function disturbs cytokinesis and results in large multinucleate or polyploid cells. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	23176
基因名	SEPTIN8

Swiss	Q92599
别名	Septin 8 (YD31950)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn