

— ABLYBIO, Help Your Research



GPA33 (YD32072) Rabbit mAb

货号: **AYD13589**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P
推荐浓度	
理论分子量	36kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse small intestine, rat small intestine
细胞定位	Membrane
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The glycoprotein encoded by this gene is a cell surface antigen that is expressed in greater than 95% of human colon cancers. The open reading frame encodes a 319-amino acid polypeptide having a putative secretory signal sequence and 3 potential glycosylation sites. The predicted mature protein has a 213-amino acid extracellular region, a single transmembrane domain, and a 62-amino acid intracellular tail. The sequence of the extracellular region contains 2 domains characteristic of the CD2 subgroup of the immunoglobulin (Ig) superfamily.
基因ID	10223
基因名	GPA33
Swiss	Q99795 (https://www.uniprot.org/uniprotkb/Q99795/entry)
别名	GPA33 (YD32072),GPA33 (YD32072) Rabbit mAb,GPA33,Glycoprotein A33

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)