

— ABLYBIO, Help Your Research



FABP7 (YD31961) Rabbit mAb

货号: **AYD13617**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	15kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG, A-549
细胞定位	Cytoplasm
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The gene encodes a small, highly conserved cytoplasmic protein that bind long-chain fatty acids and other hydrophobic ligands. The encoded protein is important in the establishment of the radial glial fiber in the developing brain. Alternative splicing and promoter usage results in multiple transcript variants encoding different isoforms. Pseudogenes of this gene are found on multiple chromosomes.
------	--

基因ID	2173
基因名	FABP7
Swiss	O15540 (https://www.uniprot.org/uniprotkb/O15540/entry)
别名	FABP7 (YD31961),FABP7 (YD31961) Rabbit mAb,FABP7,Brain lipid-binding protein,Brain-type fatty acid-binding protein,Fatty acid-binding protein 7,Mammary-derived growth inhibitor related,BLBP,FABPB,MRG

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)