

# ISG20 (YD31893) Rabbit mAb

货号: **AYD13630**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	20kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,A-549,MCF7,Mouse spleen,Rat spleen
细胞定位	Nucleus, nucleolus, Cytoplasm, Cajal body, P-body
纯化	

## 抗原信息

抗原信息	
------	--


## 靶点信息

研究背景	Interferon-induced antiviral exoribonuclease that acts on single-stranded RNA and also has minor activity towards single-stranded DNA. Exhibits antiviral activity against RNA viruses including hepatitis C virus (H CV, hepatitis A virus (HAV and yellow fever virus (YFV in an exonuclease-dependent manner. May also play additional roles in the maturation of snRNAs and rRNAs, and in ribosome biogenesis.
基因ID	3669
基因名	ISG20
Swiss	Q96AZ6

别名	ISG20 (YD31893)
----	-----------------

## 产品验证

---



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)