

— ABLYBIO, Help Your Research



Rsk 2/MAPKAP Kinase 1b (YD31679) Rabbit mAb

货号: AYD13660

产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC FC IP
推荐浓度	
理论分子量	84kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung
细胞定位	Nucleus, Cytoplasm
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 non-identical kinase catalytic domains and phosphorylates various substrates, including members of the mitogen-activated kinase (MAPK) signalling pathway. The activity of this protein has been implicated in controlling cell growth and differentiation. Mutations in this gene have been associated with Coffin-Lowry syndrome (CLS).
基因ID	6197
基因名	RPS6KA3
Swiss	P51812 (https://www.uniprot.org/uniprotkb/P51812/entry)
别名	Rsk 2/MAPKAP Kinase 1b (YD31679),Rsk 2/MAPKAP Kinase 1b (YD31679) Rabbit mAb,RPS6KA3,90 kDa ribosomal protein S6 kinase 3,Insulin-stimulated protein kinase 1,MAP kinase-activated protein kinase 1b,Ribosomal S6 kinase 2,pp90RSK2,ISP1,MAPKAPK1B,RSK2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)