

TREX1 (YD11870) Rabbit mAb

货号: **AYD13675**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	33kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HT-29
细胞定位	Nucleus, Cytoplasm, cytosol, Endoplasmic reticulum membrane
纯化	

抗原信息




抗原信息	
------	--

靶点信息

研究背景	This gene encodes a nuclear protein with 3' exonuclease activity. The encoded protein may play a role in DNA repair and serve as a proofreading function for DNA polymerase. Mutations in this gene result in Aicardi-Goutieres syndrome, chilblain lupus, Cree encephalitis, and other diseases of the immune system. Alternative splicing results in multiple transcript variants.
基因ID	11277
基因名	TREX1
Swiss	Q9NSU2

别名	TREX1 (YD11870)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn