

— ABLYBIO, Help Your Research



FOXA1 (YD31652) Rabbit mAb

货号: **AYD13676**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P ICC/IF
推荐浓度	
理论分子量	49kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat testis
细胞定位	Nucleus
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver.
------	--

基因ID	3169
基因名	FOXA1
Swiss	P55317 (https://www.uniprot.org/uniprotkb/P55317/entry)
别名	FOXA1 (YD31652),FOXA1 (YD31652) Rabbit mAb,FOXA1,Forkhead box protein A1,Transcription factor 3A,HNF3A,TCF3A

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)