

WISP1 (YD31598) Rabbit mAb

货号: **AYD13719**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	41kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	


抗原信息

抗原信息	
------	--

靶点信息

研究背景	Downstream regulator in the Wnt/Frizzled-signaling pathway (By similarity). Associated with cell survival. Adheres to skin and melanoma fibroblasts (By similarity). In vitro binding to skin fibroblasts occurs through the proteoglycans, decorin and biglycan (By similarity)
基因ID	65154
基因名	Ccn4
Swiss	Q99PP0
别名	WISP1 (YD31598)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn