

— ABLYBIO, Help Your Research



ASH2L (YD11854) Rabbit mAb

货号: **AYD13745**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC IP
推荐浓度	
理论分子量	69kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,Mouse lung,Mouse testis
细胞定位	Nucleus
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	Transcriptional regulator. Component or associated component of some histone methyltransferase complexes which regulates transcription through recruitment of those complexes to gene promoters. Component of the Set1/Ash2 histone methyltransferase (HMT complex, a complex that specifically methylates 'Lys-4' of histone H3, but not if the neighboring 'Lys-9' residue is already methylated. As part of the MLL1/MLL complex it is involved in methylation and dimethylation at 'Lys-4' of histone H3. May play a role in hematopoiesis. In association with RBBP5 and WDR5, stimulates the histone methyltransferase activities of KMT2A, KMT2B, KMT2C, KMT2D, SETD1A and SETD1B.
基因ID	9070
基因名	ASH2L
Swiss	Q9UBL3 (https://www.uniprot.org/uniprotkb/Q9UBL3/entry)
别名	ASH2L (YD11854),ASH2L (YD11854) Rabbit mAb,ASH2L,ASH2-like protein,ASH2L1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)