

Annexin A1/ANXA1 (YD31173) Rabbit mAb

货号: **AYD13840**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	39kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C2C12,C6
细胞定位	Nucleus, Cytoplasm, Cell projection, cilium, Cell membrane, Membrane, Endosome membrane, Basolateral cell membrane, Apical cell membrane, Lateral cell membrane, Secreted, extracellular space, extracellular exosome, Cytoplasmic vesicle, secretory vesicle lumen, phagocytic cup, Early endosome, Cytoplasmic vesicle membrane
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a membrane-localized protein that binds phospholipids. This protein inhibits phospholipase A2 and has anti-inflammatory activity. Loss of function or expression of this gene has been detected in multiple tumors.
基因ID	301
基因名	ANXA1

Swiss	P04083
别名	Annexin A1/ANXA1 (YD31173)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn