

# Phospho-TrkA (Tyr490)/TrkB (Tyr516)/TrkC (Tyr516) (YD11627) Rabbit mAb

货号: AYD13864

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	92kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell membrane, Endosome membrane, Early endosome membrane, Cell projection, axon, dendrite, Cytoplasm, perinuclear region, Postsynaptic density
纯化	

## 抗原信息

抗原信息	
------	--


## 靶点信息

研究背景	This gene encodes a member of the neurotrophic tyrosine receptor kinase (NTRK) family. This kinase is a membrane-bound receptor that, upon neurotrophin binding, phosphorylates itself and members of the MAPK pathway. Signalling through this kinase leads to cell differentiation. Mutations in this gene have been associated with obesity and mood disorders. Alternative splicing results in multiple transcript variants.
基因ID	4915
基因名	NTRK2

<b>Swiss</b>	Q16620
别名	Phospho-TrkA (Tyr490)/TrkB (Tyr516)/TrkC (Tyr516) (YD11627)

## 产品验证

---



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)