

# Phospho-PSD95 (Ser418) (YD21006) Rabbit mAb

货号: AYD13885

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	80kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Mouse brain,Rat brain
细胞定位	Cell membrane, Postsynaptic density, Synapse, Cytoplasm, Cell projection, axon, dendritic spine, dendrite, Presynapse
纯化	

## 抗原信息

抗原信息	
------	--

## 靶点信息


研究背景	This gene encodes a member of the membrane-associated guanylate kinase (MAGUK) family. It heteromultimerizes with another MAGUK protein, DLG2, and is recruited into NMDA receptor and potassium channel clusters. These two MAGUK proteins may interact at postsynaptic sites to form a multimeric scaffold for the clustering of receptors, ion channels, and associated signaling proteins. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	1742

基因名	DLG4
Swiss	P78352
别名	Phospho-PSD95 (Ser418) (YD21006)

## 产品验证

---

Mouse brain



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)