

— ABLYBIO, Help Your Research



Phospho-Bcl-2 (Ser70) (YD20039) Rabbit mAb

货号: **AYD13900**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC IP
推荐浓度	
理论分子量	26kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Mitochondrion outer membrane, Nucleus membrane, Endoplasmic reticulum membrane, Cytoplasm
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants.
------	---

基因ID	596
基因名	BCL2
Swiss	P10415 (https://www.uniprot.org/uniprotkb/P10415/entry)
别名	Phospho-Bcl-2 (Ser70) (YD20039),Phospho-Bcl-2 (Ser70) (YD20039) Rabbit mAb,BCL2

产品验证

Mouse brain

kDa

250 -
150 -
100 -
70 -
50 -
40 -
35 -
25 -
20 -
15 -

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)