

— ABLYBIO, Help Your Research



Phospho-PDGFR alpha (Tyr720) (YD19131) Rabbit mAb

货号: AYD13909

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	123kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C6,Mouse lung
细胞定位	Cell membrane, Cell projection, cilium, Golgi apparatus
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a cell surface tyrosine kinase receptor for members of the platelet-derived growth factor family. These growth factors are mitogens for cells of mesenchymal origin. The identity of the growth factor bound to a receptor monomer determines whether the functional receptor is a homodimer or a heterodimer, composed of both platelet-derived growth factor receptor alpha and beta polypeptides. Studies suggest that this gene plays a role in organ development, wound healing, and tumor progression. Mutations in this gene have been associated with idiopathic hypereosinophilic syndrome, somatic and familial gastrointestinal stromal tumors, and a variety of other cancers.
基因ID	5156
基因名	PDGFRA
Swiss	P16234 (https://www.uniprot.org/uniprotkb/P16234/entry)
别名	Phospho-PDGFR alpha (Tyr720) (YD19131), Phospho-PDGFR alpha (Tyr720) (YD19131) Rabbit mAb, PDGFR A, Alpha platelet-derived growth factor receptor, Alpha-type platelet-derived growth factor receptor, CD140 antigen-like family member A, CD140a antigen, Platelet-derived growth factor alpha receptor

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)