

— ABLYBIO, Help Your Research



FXR1 (YD31258) Rabbit mAb

货号: **AYD13931**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	70kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3 cells, HEK-293 cells, HeLa cells, mouse heart tissue
细胞定位	Cytoplasm, Cytoplasmic ribonucleoprotein granule, Stress granule, Cell projection, dendrite, dendritic spine, axon, Nucleus envelope, Postsynapse
纯化	亲和纯化




抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is an RNA binding protein that interacts with the functionally-similar proteins FMR1 and FXR2. These proteins shuttle between the nucleus and cytoplasm and associate with poly ribosomes, predominantly with the 60S ribosomal subunit. Three transcript variants encoding different isoforms have been found for this gene.
基因ID	8087
基因名	FXR1
Swiss	P51114 (https://www.uniprot.org/uniprotkb/P51114/entry)
别名	FXR1 (YD31258),FXR1 (YD31258) Rabbit mAb,FXR1,FMR1 autosomal homolog 1,hFXR1p

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)