

— ABLYBIO, Help Your Research



# Phospho-CD3 zeta (Tyr83) (YD16224) Rabbit mAb

货号: AYD13935

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF FC IP
推荐浓度	
理论分子量	19kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse thymus
细胞定位	Cell membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is T-cell receptor zeta, which together with T-cell receptor alpha/beta and gamma/delta heterodimers, and with CD3-gamma, -delta and -epsilon, forms the T-cell receptor-CD3 complex. The zeta chain plays an important role in coupling antigen recognition to several intracellular signal-transduction pathways. Low expression of the antigen results in impaired immune response. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
基因ID	919
基因名	CD247
Swiss	P20963 ( <a href="https://www.uniprot.org/uniprotkb/P20963/entry">https://www.uniprot.org/uniprotkb/P20963/entry</a> )
别名	Phospho-CD3 zeta (Tyr83) (YD16224),Phospho-CD3 zeta (Tyr83) (YD16224) Rabbit mAb,CD247,T-cell receptor T3 zeta chain,CD3Z,T3Z,TCRZ

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)