

— ABLYBIO, Help Your Research



Phospho-c-Myc (Thr58) (YD16033) Rabbit mAb

货号: **AYD13938**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB ICC/IF FC |
| 推荐浓度 | |
| 理论分子量 | 51kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | 293T |
| 细胞定位 | Nucleus, nucleoplasm, nucleolus, Cytoplasm, Chromosome |
| 纯化 | 亲和纯化 |

抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|---|
| 研究背景 | The protein encoded by this gene is a multifunctional, nuclear phosphoprotein that plays a role in cell cycle progression, apoptosis and cellular transformation. It functions as a transcription factor that regulates transcription of specific target genes. Mutations, overexpression, rearrangement and translocation of this gene have been associated with a variety of hematopoietic tumors, leukemias and lymphomas, including Burkitt lymphoma. There is evidence to show that alternative translation initiations from an upstream, in-frame non-AUG (CUG) and a downstream AUG start site result in the production of two isoforms with distinct N-termini. The synthesis of non-AUG initiated protein is suppressed in Burkitt's lymphomas, suggesting its importance in the normal function of this gene. |
| 基因ID | 4609 |
| 基因名 | MYC |
| Swiss | P01106 (https://www.uniprot.org/uniprotkb/P01106/entry) |
| 别名 | Phospho-c-Myc (Thr58) (YD16033),Phospho-c-Myc (Thr58) (YD16033) Rabbit mAb,MYC,Class E basic helix-loop-helix protein 39,Proto-oncogene c-Myc,Transcription factor p64,BHLHE39 |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)