

Thioredoxin/TRX (YD11891) Rabbit mAb

货号: **AYD13952**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	12kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,H460,MCF7,HT-29,SKOV3
细胞定位	Nucleus, Cytoplasm, Secreted
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene acts as a homodimer and is involved in many redox reactions. The encoded protein is active in the reversible S-nitrosylation of cysteines in certain proteins, which is part of the response to intracellular nitric oxide. This protein is found in the cytoplasm. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	7295
基因名	TXN
Swiss	P10599

别名	Thioredoxin/TRX (YD11891)
----	---------------------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn