

— ABLYBIO, Help Your Research



P70 S6 Kinase beta (YD14930) Rabbit mAb

货号: **AYD14017**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	53kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SH-SY5Y,HL-60,SW620,MCF7,Mouse liver,Mouse kidney,Mouse spleen,Rat liver
细胞定位	Cytoplasm, Nucleus
纯化	亲和纯化

抗原信息





抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains a kinase catalytic domain and phosphorylates the S6 ribosomal protein and eukaryotic translation initiation factor 4B (eIF4B). Phosphorylation of S6 leads to an increase in protein synthesis and cell proliferation.
------	---

基因ID	6199
基因名	RPS6KB2
Swiss	Q9UBS0 (https://www.uniprot.org/uniprotkb/Q9UBS0/entry)
别名	P70 S6 Kinase beta (YD14930),P70 S6 Kinase beta (YD14930) Rabbit mAb,RPS6KB2,70 kDa ribosomal protein S6 kinase 2,S6 kinase-related kinase,Serine/threonine-protein kinase 14B,p70 ribosomal S6 kinase beta,STK14B

产品验证



实验步骤

