

— ABLYBIO, Help Your Research



Phospho-Filamin A (Ser1083) (YD13051) Rabbit mAb

货号: **AYD14041**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	281kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,U-251MG,Mouse lung,Mouse spleen
细胞定位	Cytoplasm, cell cortex, cytoskeleton, Perikaryon, Cell projection, growth cone, podosome
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is an actin-binding protein that crosslinks actin filaments and links actin filaments to membrane glycoproteins. The encoded protein is involved in remodeling the cytoskeleton to effect changes in cell shape and migration. This protein interacts with integrins, transmembrane receptor complexes, and second messengers. Defects in this gene are a cause of several syndromes, including periventricular nodular heterotopias (PVNH1, PVNH4), otopalatodigital syndromes (OPD1, OPD2), frontometaphyseal dysplasia (FMD), Melnick-Needles syndrome (MNS), and X-linked congenital idiopathic intestinal pseudoobstruction (CIIPX). Two transcript variants encoding different isoforms have been found for this gene.
基因ID	2316
基因名	FLNA
Swiss	P21333 (https://www.uniprot.org/uniprotkb/P21333/entry)
别名	Phospho-Filamin A (Ser1083) (YD13051), Phospho-Filamin A (Ser1083) (YD13051) Rabbit mAb, FLNA, Actin-binding protein 280, Alpha-filamin, Endothelial actin-binding protein, Filamin-1, Non-muscle filamin, FLN, FLN1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)