

Phospho-delta 1 Catenin (Tyr96) (YD12918) Rabbit mAb

货号: **AYD14051**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	108kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,NIH/3T3,C6
细胞定位	Cell junction, adherens junction, Cytoplasm, Nucleus, Cell membrane
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length nature of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thiodoxin-related transmembrane protein 2 (TMX2) gene.
基因ID	1500
基因名	CTNND1

Swiss	O60716
别名	Phospho-delta 1 Catenin (Tyr96) (YD12918)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn